



PolyStain DS Kit - for 2 Mouse antibody on Human tissue

For co-localization (Emerald/Permanent Red)

NB-23-00099- 3(120 ml)

NB-23-00099- 2(36 ml)

NB-23-00099- 1(12 ml)

PolyStain DS Kit - for 2 Mouse antibody on Human tissue For co-localization (Emerald/Permanent Red)

NB-23-00099-1; NB-23-00099-2; NB-23-00099-3

Storage: 2-8°C

INTENDED USE:

PolyStain DS Kit is designed to use with user supplied two mouse primary antibodies to detect two distinct antigens on human tissue or cell samples. This kit has been tested in paraffin tissue. However, this kit can be used on frozen specimen and freshly prepared monolayer cell smears. Double staining is one of the most common methods used in immunohistochemistry to evaluate two distinct antigens in a single tissue.

PolyStain DS Kit from NeoBiotech Labs supplies two polymer enzyme conjugates: HRP polymer anti-Mouse IgG and AP polymer anti-Mouse IgG with two distinct substrates/chromogens, Emerald (Green color, use with HRP polymer anti-Mouse IgG) and Permanent-Red (Red color, use with AP polymer anti-Mouse IgG).

Simplified steps offer a much faster protocol and the blocking buffers prevent false positives when using two primary antibodies from the same host species. Another advantage of NeoBiotech Kit, it allows the researcher to visualize when two proteins are co-localized because of the color change when the chromogens overlap that can be semi-quantitative. For example, if the area of co-localization stains blue, the antigen indicated by Emerald is at higher concentration in the cell and if the color is purple than the antigen indicated by Permanent-Red is expressed at higher concentrations. PolyStain DS Kit is non-biotin system that avoids endogenous biotin nonspecific binding.

KIT COMPONENTS:

Component No.	Content	12mL Kit	36mL Kit	120mL Kit
Reagent 1	AP Polymer anti-Mouse IgG (RTU)	6ml	18ml	60ml
Reagent 2A	Permanent Red Substrate (RTU)	7ml	18ml	60ml
Reagent 2B	Permanent Red Activator (5x)	1.4ml	3.6ml	12ml
Reagent 2C	Permanent Red Chromogen (100x)	70µL	180µL	0.6ml
Reagent 3	Antibody Blocker (40x)	2x15mL	50mL	125mL
Reagent 4A	DS-MM Blocker A (RTU)	6mL	18mL	60mL
Reagent 4B	DS-MM Blocker B (RTU)	6mL	18mL	60mL
Reagent 5	HRP Polymer anti-Mouse IgG (RTU)	6mL	18mL	60mL
Reagent 6	Emerald Chromogen (RTU)	7mL	18mL	70mL
Reagent 7	U-Mount (RTU)	3mL	9mL	NA

RECOMMENDED PROTOCOL:

1. Fixation: To ensure the quality of the staining and obtain reproducible performance, user needs to supply appropriately fixed tissue and well prepared slides.
2. Tissue need to be adhered to the slide tightly to avoid tissue falling off.
3. Paraffin embedded section must be deparaffinized with xylene and rehydrated with a graded series of ethanol before staining.
4. Cell smear samples should be made as much monolayer as possible to obtain satisfactory results.
5. Three control slides will aid the interpretation of the result: positive tissue control, reagent control (slides treated with Isotype control reagent), and negative control.
6. Proceed IHC staining: **DO NOT** let specimen or tissue dry from this point on.

Reagent	Staining Procedure	Incubation Time (Min.)
1. Peroxidase and alkaline phosphatase Blocking Reagent Not provided	a. Incubate slides in peroxidase and alkaline phosphatase blocking reagent. We recommend NeoPure Dual Enzyme Block NB-23-00193 for 10 minutes... b. Rinse the slide using distilled water.	10 min.
2. HIER Pretreatment: Refer to antibody data sheet.	a. Heat Induced Epitope Retrieval (HIER) may be required for primary antibody suggested by vendor. b. Wash with PBS containing 0.05% Tween-20 for 2 min., 3 times.	
3. Preblock (optional)	For paraffin section, Improved polymer formula saves the need for a pre-block step. However some primary antibodies may still require it, this information should be determined by the user prior to using DS kit. For paraffin section, Improved formula saves the need for a preblock step. For frozen tissue, preblock may or may not be required depending on fixative. (Preblock catalogue No. NB-23-00169 was Recommended.)	
4. Mouse Antibody 1: Supplied by user	Notes: Investigator needs to optimize dilution and incubation times prior to double staining. a. Apply 2 drops or enough volume of mouse primary antibody 1 to cover the tissue completely. Incubate in moist chamber for 30-60 min. b. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times	30 - 60 min.

<p>5. Reagent 1 AP polymer anti-Mouse IgG (RTU)</p>	<p>a. Apply 2 drops (100µL) of Reagent 1 (AP polymer anti-Mouse IgG) to cover each section. b. Incubate in moist chamber for 15 min. c. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times. d. Rinse well with distilled or tap water.</p>	<p>15 min</p>
<p>6. Reagent 2A, 2B, 2C Reagent 2A: Permanent Red Substrate (RTU) Reagent 2B: Permanent Red Activator (5x) Reagent 2C: Permanent Red Chromogen (100x)</p>	<p>a. Add 200µL of Reagent 2B (Activator) into 1mL of Reagent 2A (Substrate buffer) and mix well. Add 10µL of Reagent 6C (Chromogen) into the mixture and mix well. [Note: For fewer slides, Add 100µL of Reagent 2B (Activator) into 500µL of Reagent 2A (Substrate buffer) and mix well. Add 5µL of Reagent 2C (Chromogen) into the mixture and mix well.] b. Apply 2 drops (100µL) or enough volume of Permanent Red working solution to completely cover the tissue. Incubate for 10 min, observe appropriate color development. c. Rinse well with distilled water.</p>	<p>10 min</p>
<p>7. Reagent 3: Antibody Blocker(40x) (Optional) Must test if antibody/antigen interaction is heat sensitive. Please skip this step if antigen retrieval is used for 2nd Ms Primary Antibody.</p>	<p>Note: This step will block antibodies of previous step so no cross reaction will occur at end of protocol. HIER can be done immediately after Antibody Blocker step if only one primary antibody requires antigen retrieval. For frozen tissue a lower temperature of 65°C must be used for Antibody Blocker (Reagent 3) to prevent tissue from dissociating from slide. a. Use hot plate or water bath to heat diluted Reagent 3 to 1x solution (1 part of Antibody Blocker in 39 parts of distilled water) to 80°C. Make enough volume to cover the tissue in beaker. b. Put slides in heated Antibody Blocker for 10 minutes at 80°C. c. Remove slides from the Antibody blocker; cool slides 5 seconds. d. Rinse slides in multiple changes of distilled water. e. Wash with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</p>	<p>10 min</p>
<p>8. Reagent 4A DS-MM Blocker A</p>	<p>a. Apply 2 drops or enough volume of Reagent 4A (DSMM Blocker A) to cover the tissue completely. Mix well on the slide and Incubate in moist chamber for 30 min. b. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times</p>	<p>30 min</p>
<p>9. Reagent 4B DS-MM Blocker B</p>	<p>a. Apply 2 drops or enough volume of Reagent 4B (DSMM Blocker B) to cover the tissue completely. Mix well on the slide and Incubate in moist chamber for 10 min. b. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times.</p>	<p>5 min.</p>
<p>10. Mouse antibody 2: Supplied by user</p>	<p>Notes: Investigator needs to optimize dilution and incubation times prior to double staining. a. Apply 2 drops or enough volume of mouse primary antibody 2 to cover the tissue completely. b. Rinse with PBS containing 0.05% Tween-20 for 2 min., 3 times</p>	<p>30- 60 min.</p>

<p>11. Reagent 5 HRP polymer anti-Mouse IgG (RTU)</p>	<ol style="list-style-type: none"> Apply 2 drops (100µL) of Reagent 5 AP polymer AntiMouse IgG to cover each section. Incubate in moist chamber for 15 min. Rinse with PBS containing 0.05% Tween-20 for 2 min, 3 times. Rinse with tap water. 	<p>15 min</p>
<p>12. Counterstain (Optional)</p> <p>Not provided</p>	<ol style="list-style-type: none"> Dip the slide in diluted hematoxylin for 5 seconds. (You may dilute hematoxylin 1:5 in dH₂O). DO NOT over stain with hematoxylin. Rinse thoroughly with tap water for 2min. Put slides in PBS for 5 seconds to blue, DO NOT over blue. Rinse well in distilled or tap water for 2min. Wash with PBS/0.05% Tween20 for 2 min, 3 times. 	<p>5 seconds</p>
<p>13. Reagent 6 Emerald Chromogen(RTU)</p>	<ol style="list-style-type: none"> Apply 1 to 2 drops (50-100µL) of Reagent 6 (Emerald Chromogen) to cover the tissue completely. Incubate in moist chamber for 5 minutes. Wash slides in tap water for 1 minute. Rinse with distilled water. <p>Important to READ: <i>Emerald Chromogen is water soluble, do counter stain first. Do not leave slides sitting in water. Always stain Emerald chromogen AFTER Permanent Red stain because GBI-Permanent Red removes the Emerald and after hematoxylin.</i></p>	<p>5 min</p>
<p>14. Dehydrate section</p>	<p>Note: Please wipe off extra water and air dry slides before dehydration and clear.</p> <ol style="list-style-type: none"> Dehydrate with 85% ethanol 20seconds. Dehydrate with 95% ethanol 20seconds. Dehydrate with 100% ethanol 20seconds. Dehydrate with 100% ethanol 20seconds. Dehydrate with 100% ethanol 20seconds. Dehydrate with xylene 20seconds. <p>CAUTION: DO NOT dehydrate with xylene longer than 20 seconds! It will erase Permanent Red stain!</p>	
<p>15. Reagent 7 U-Mount(RTU)</p>	<ol style="list-style-type: none"> Apply 1 to 2 drops (50-100µL) of Reagent 7 (U Mount) to cover the tissue section and apply glass coverslip. Apply force to coverslip to squeeze out any extra mountant and bubbles for optimal clarity. Removing excess also to prevent leaching of Permanent Red stain. 	

TROUBLE SHOOT:

PROBLEM	TIPS
Uneven stain on 2 primary antibodies	<ol style="list-style-type: none"> 1. Need to adjust the titer of each antibody. 2. The amount of each protein expressed on tissue may be different. 3. Set slides in water too long so that Emerald is washed away. 4. Set slides in Xylene too long so that Permanent Red is washed away
Emerald Chromogen is blue not green when non co-localized with Permanent Red.	Emerald should be green when not co-localized with Permanent Red. If Emerald chromogen is blue the titer on the primary antibody is not dilute enough for the protocol. Re-titer primary antibodies individually first.
No stain on 1 or 2 antibodies	Missing steps or step reversed.
Green Background on the slide	<ol style="list-style-type: none"> 1. Titer primary antibody. 2. Use 10% Donkey serum, goat or horse serum as a preblock.
Permanent Red is leaching	<ol style="list-style-type: none"> 1. Use fresh 100% ethanol and xylene. 2. Slide sat too long in xylene. Do not go over 20seconds!
Artifacts on slides	Slides not completely dried before mount. Use fresh 100% Ethanol and xylene.

PRECAUTIONS:

DAB may be carcinogenic. Please wear gloves and take other necessary precautions.

FOR RESEARCH USE ONLY

Work Sheet for NB-23-00099 Kit

We designed this work sheet to help you track of each step. We recommend you use this sheet to record the actual time of each step conducted as it will be helpful for questions with our technical support. To insure that all steps are done properly, we recommend that the user fill in the actual time of their experimental step and any variation. Results will vary if time recommendations are not followed.

RTU translates to ready to use.

- Used for tester to check “√ “each step during the experiment
- Steps follow after de-paraffinization
- Refer to insert for details of each step

NB-23-00099 Protocol is suitable when both mouse and rabbit primary antibodies need or do not need pre-treatment step

Protocol Step	NB-23-00097 Protocol	Experiment 1 Date:	Experiment 2 Date:	Experiment 3 Date:	Experiment 4 Date:
Step 1	Peroxidase & Alkaline Phosphatase Block User supplied				
Step 2 (Optional)	HIER if needed User supplied (up to 60 min)				
Step 3 Optional	Preblock if needed User supplied				
Step 4	Mouse 1°Ab #1 User supplied (30-60 min.)				
Step 5	Reagent 1 AP Polymer anti-Mouse IgG RTU(15min) Rinse with distilled water				
Step 6	Reagent 2A, Reagent 2B & Reagent 2C Permanent Red requires mixing (10min)				

Step 7	Reagent 3 Antibody Blocker (10 min)				
Step 8	Reagent 4A DS-MM Blocker A RTU (30 min)				
Step 9	Reagent 4B DS-MM Blocker B RTU (5 min)				
Step 10	Ms 1°Ab #2 User supplied (30-60 min)				
Step 11	Reagent 5 HRP Polymer anti-Mouse IgG RTU (15min) Rinse with tap water				
Step 12	Counter stain (Do not over counter stain) Hematoxylin User supply Wash with PBS/0.05% Tween20 for 2 min, 3 times				
Step 13	Reagent 6 Emerald Chromogen RTU (5min)				
Step 14	Dehydrate section 20seconds for each step It is important to follow the protocol.				
Step 15	Reagent 7 U-Mount RTU Mount & coverslip				
Result	Stain pattern on controls are correct: Fill in Yes or NO				