

## PCR Enhancer

**Ref: NB-03-0153 500 µl**

### Description

PCR Enhancer System is used in conjunction with DNA polymerases to optimize PCR of problematic and/or GC-rich templates. The system offers higher primer specificity, broader magnesium concentration, broader annealing temperature, and improved thermostability DNA polymerases.

### Usage

Add 2 – 8 µl PCR Enhancer to 25 µl amplification system, the other amplification conditions and procedure remain the same.

**Store at -20°C**

### PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.