

## Monoclonal anti RFP / DsRed (*Red Fluorescent Protein*)

#Cat: NB-01-0107

Size: 100 µl

**Clones:** 1DS-1A6

**Immunogen:** RFP recombinant protein

**Murine Isotype:** IgG1, k

**Specificity:** RFP, RFP-fusion proteins

**Applications:**

- ELISA
- Western Blot
- Immunocytochemistry

**Recommended dilutions:** 1/500-1/5000

**Storage:**

- Store at -20°C (up to 3 years)
- Avoid repeat freezing and thawing cycles. Does not contain any preservative

**These products are intended for research purposes only.**

**MATERIAL SAFETY DATA SHEET INFORMATION**

To the best of our knowledge, these products do not require a Material Safety Data Sheet. However, all the properties of these products (and, if applicable, each of their components) have not been thoroughly investigated. Therefore, we recommend that you use gloves and eye protection, and wear a laboratory coat when working with these products.