

## Rabbit anti-Human Neurotensin Receptor 1 (NTR1) (AA 398-418) (CT)

#Cat: NB-49-005-50ul

Size: 50ul

<b>Product Description:</b>	NT1 antibody, human, C-term. NTR1 Neurotensin 1 receptor G Protein-coupled receptor Polyclonal antibody raised against the C-terminus of the human NT1 receptor
<b>Immunogen:</b>	residues 398-418 coupled to KLH (NP_002522) KLH- kadsvssnhtlssnatretly
<b>Host:</b>	Rabbit
<b>Form:</b>	serum with 0.02% NaN <sub>3</sub> as preservative
<b>Quantity:</b>	50 µl
<b>Applications:</b>	IC, W W = Western blot (1:5000) IC = Immunocytochemistry (1:5000)
<b>Species cross-reactivity:</b>	human
<b>Sensitivity:</b>	50 ng
<b>Storage:</b>	4 °C short term -20 °C long term

