

Rabbit anti-Human Bombesin-Like Receptor 3 (BRS-3) (AA 342-361) (CT)

#Cat: NB-49-010-50ul

Size: 50ul

Product Description:	BRS-3 antibody, human, C-term. Bombesin-like receptor 3 Bombesin receptor G Protein-coupled receptor Polyclonal antibody raised against the C-terminus of the human BRS-3 receptor
Immunogen:	residues 342-361coupled to KLH (NP_001718) KLH-seisvtsftgcsvkqaedrf
Host:	Rabbit
Form:	serum with 0.02% NaN3 as preservative
Quantity:	50 µl
Applications:	ICH (p), IC, W W = Western blot (1:5000) IC = Immunocytochemistry (1:5000) ICH (p) = Immunohistochemistry of formalin-fixed, paraffin- embedded tissues (1:2500)
Species cross-reactivity:	human
Sensitivity:	50 ng
Storage:	4° C short term -20° C long term



BRS-3 immunoreactivity in human pancreatic carcinoma.