

Nutrient Agar

(ISO 6579, ISO 10273, ISO 19250)

(Dehydrated Culture Media) for microbiology

Cat # NB-42-74553-500g Size: 500g

Composition (g/l)

Meat Extract	3,0
Gelatine Peptone	5,0
Agar	15,0

pH: 6.8 +/- 0.2

Physical Data:

Dehydrated culture media of light beige colour.

Notes

Medium for the differential subculture of *Pseudomonas aeruginosa* according to UNE-EN 1278:22 of water.

Preparation

Make a suspension of 23 g in 1 L of distilled water. Heat stirring constantly to boiling point and until completely dissolved. Boil for 1 minute. Sterilize at 121°C for 15 minutes.

Storage

Storage in a cool and dry place. Once opened keep powdered medium closed to avoid hydration

For reference only

For Research Use Only. Not for Diagnostic or Therapeutic Use.