
Rabbit anti-Human Somatostatin Receptor SST1 (SSTR1)

#Cat: NB-49-019-50ul

Size: 50ul

Batch:	S15-11-I
Quantity:	50 µl
Host:	Rabbit
Immunogen:	Peptide of human SST1 Receptor coupled to KLH
Presentation:	Whole Antiserum, liquid
Specificity:	This antiserum recognizes the C-terminus of human SST1 receptor. Specificity of antibody binding was accomplished by blocking with is cognate peptide.
Application:	Immunolocalization
Titer:	Working dilution has been determined to be 1:3000 – 1:4000 in Immunohistochemistry.
Storage:	Aliquot and store at -20°C or colder. Avoid repeated freeze and thawing cycles.