
Rabbit anti-Human Somatostatin Receptor SST4 (SSTR4)

#Cat: NB-49-024-50ul

Size: 50ul

Batch: S15-11-I**Quantity:** 50 µl**Host:** Rabbit**Immunogen:** Peptide of human SST4 Receptor coupled to KLH**Presentation:** Whole Antiserum, liquid**Specificity:** This antiserum recognizes the C-terminus of human SST4 receptor. Specificity of antibody binding was accomplished by blocking with its cognate peptide.**Application:** Immunolocalization of SST4 receptor**Titer:** Working dilution has been determined to be 1:3000 – 1:4000 in Immunohistochemistry.**Storage:** Aliquot and store at -20°C or colder. Avoid repeated freeze and thawing cycles.