
Rabbit anti-Mouse Somatostatin Receptor SST3 (SSTR3)

#Cat: NB-49-017-50ul

Size: 50ul

Batch:	S15-11-I
Quantity:	50 µl
Host:	Rabbit
Immunogen:	Peptide of rat and mouse SSTR3 Receptor coupled to KLH
Presentation:	Whole Antiserum, liquid
Reconstitution:	Reconstitute in an appropriate volume of distilled water and store in aliquots at -20°C or colder.
Specificity:	This antiserum recognizes the C-terminus of rat and mouse SSTR3 receptor. Specificity of antibody binding was accomplished by blocking with its cognate peptide.
Application:	Immunolocalization
Titer:	Working dilution has been determined to be 1:3000 – 1:4000 in Immunohistochemistry.
Storage:	Aliquot and store at -20°C or colder. Avoid repeated freeze and thawing cycles.