

# Dry Bath Incubator - N70 Series

The N70 dry block heating is designed for testing milk for antibacterial substances such as antibiotics and sulphonamide residues. It includes a digital thermometer and timer. Heater is either AC or DC compatible. AC adapter is included, DC adapter for vehicle must be purchased separately.

## Features :

- A personal dry block heating and for digital dry block heating
- It pretty much fits anywhere on your lab bench
- It has a very simple touch pad digital display for temperature and time selection
- It's only 4-inch wide
- High quality -- low price



## Specifications :

Model	N70
Cat#	NB-12-0011
Temperature Control Range	RT+5°C ~ 80°C
Stability at 80°C	±0.5°C
Uniformity at 40°C	±0.3°C
Display Accuracy / Display	0.1°C / LCD
Heating Time (25 to 80°C)	≤12 min
Time Range	1min ~ 99h59min
Dimension LxWxH (mm)	135 x 99 x 33,5
Weight (kg)	3.0