

## Oligo d(T)<sub>15</sub> Primer

**Ref: NB-03-0151 20 µl**

### Description

Oligo d(T)<sub>15</sub> Primer This primer is a synthetic single-stranded 15-mer oligonucleotide with 5'- and 3'-hydroxyl ends. This primer is supplied as an aqueous solution and can be used for first strand cDNA synthesis.

### Application

First strand cDNA synthesis.

### Concentration

50 µM

### Quality Control

Functionally tested in first strand cDNA synthesis.

### Storage

Store at -20°C.