

Human Neuroligin-4 X linked protein (Recombinant)

#Cat : NB-22-77838-20	Size: 20µg
#Cat : NB-22-77838-100	Size: 100µg
#Cat : NB-22-77838-500	Size: 500µg

General Information

Product Type	Proteins
Short Description	Recombinant-Human Neuroligin-4 X linked-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

Product Properties

Concentration	
Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Dilution Range	>90%, as determined by SDS-PAGE and HPLC NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

Target Information

Gene ID	57502
Gene Symbol	NLGN4X
Uniprot ID	NLGNX_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.