

EPI (Epinephrine/Adrenaline) Competitive ELISA Kit

#Cat: NB-22-46307-24	Size: 24T
#Cat: NB-22-46307-48	Size: 48T
#Cat: NB-22-46307-96	Size: 96T

General Information

Product Type	Kits
Short Description	This EPI Competitive ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in all serum, plasma and other biological fluids.
Application	ELISA
Reactivity	All

Product Properties


Detection	Colormetric
Concentration	18.75pg/mL
Detection Limit	31.25~2000pg/mL
Assay Time	2.5h

Storage Instructions If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

Target Information

Sample Types	Serum, plasma and other biological fluids.
Gene Symbol	
Uniprot ID	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence.

This product must not be used as for diagnostic or other medical purposes.