

Neo DMEM High Glucose (4.5 g/l), w/o L-Glutamine, w/o Phenol Red
#Cat: NB-58-0020 Size: 500ml
Product Specifications

Appearance	Clear solution
CO₂ concentration, optimum	8.5 %
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping conditions	Ambient

Formulation

Components	Concentration mg/L
Amino Acids:	
Glycine	30.00
L-Arginine HCl	84.00
L-Cystine 2 HCl	62.57
L-Histidine HCl H ₂ O	42.00
L-Isoleucine	105.00
L-Leucine	105.00
L-Lysine HCl	146.00
L-Methionine	30.00
L-Phenylalanine	66.00
L-Serine	42.00
L-Threonine	95.00
L-Tryptophan	16.00
L-Tyrosine 2 Na 2 H ₂ O	103.79
L-Valine	94.00
Vitamins:	
Choline chloride	4.00
D-Calcium Pantothenate	4.00
Folic Acid	4.00
myo-Inositol	7.20
Nicotinamide	4.00
Pyridoxal HCl	4.00
Riboflavin	0.40
Thiamine HCl	4.00

Inorganic Salts:	
CaCl ₂ 2 H ₂ O	265.00
Fe (NO ₃) ₃ 9 H ₂ O	0.10
KCl	400.00
MgSO ₄ 7H ₂ O	200.00
NaCl	6400.00
NaHCO ₃	3700.00
NaH ₂ PO ₄ 2H ₂ O	141.73
Other Components:	
D-Glucose	4500.00

Precautions and Disclaimer

This product is for research use only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (info@neo-biotech.com) or phone (+33 9 77 40 09 09).