

Human IL6 (Interleukin-6) Sandwich ELISA Kit

#Cat: NB-22-71813-24	Size: 24T
#Cat: NB-22-71813-48	Size: 48T
#Cat: NB-22-71813-96	Size: 96T

General Information


Product Type	Kits
Short Description	This IL6 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids
Application	ELISA
Reactivity	Human

Product Properties

Detection	Colormetric
Sensitivity	0.94pg/mL
Detection Limit	1.56~100pg/mL
Assay Time	3.5h
Storage Instructions	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

Target Information

Sample Types	serum, plasma and other biological fluids
Gene Symbol	IL6
Uniprot ID	IL6_HUMAN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence.
This product must not be used as for diagnostic or other medical purposes.