

Rabbit anti-Human Growth Hormone Releasing Hormone Receptor (AA 403-422) (CT)

#Cat: NB-49-011-50ul

Size: 50ul

Product Description:	GHRH receptor antibody, human, C-term. Growth hormone releasing hormone receptor G Protein-coupled receptor Polyclonal antibody raised against the C-terminus of the human GHRH receptor
Immunogen:	residues 403-422 coupled to KLH (NP_000814) KLH- RTRAKWTTSPRSAAKVLTSM
Host:	Rabbit
Form:	serum with 0.02% NaN ₃ as preservative
Quantity:	50 µl
Applications:	ICH (p), IC, W W = Western blot (1:5000) IC = Immunocytochemistry (1:5000) ICH (p) = Immunohistochemistry of formalin-fixed, paraffin- embedded tissues (1:2500)
Species cross-reactivity:	human
Sensitivity:	50 ng
Storage:	4 ⁰ C short term -20 ⁰ C long term



GHRH receptor immunoreactivity in human ovarian carcinoma.