


Rabbit anti-Human Somatostatin Receptor SST2A (SSTR2A)

Antibody (affinity-purified)

#Cat: NB-49-015-1ml Size: 1ml

#Cat: NB-49-015-500ul Size: 500µl



SSTR2A immunoreactivity in a human neuroendocrine tumor.

Description

Rabbit anti-Somatostatin Receptor 2A (SSTR2A), Human, Rat, Mouse

Product Information

Host:	Rabbit
Applications:	IHC (p), IC, WB, IP
Reactivity:	Human, Rat, Mouse
Clonality:	Monoclonal
Conjugation:	Unconjugated
Isotype:	IgG
Purification:	affinity purified
Dilution Range:	Western blot (1:1000) Immunocytochemistry (1:500) Immunohistochemistry of formalin-fixed, paraffinembedded tissues (1:100)
Storage Instruction:	4°C (do not freeze)

Target

Immunogen:	Amino acid sequence 355-369 (ETQRTLLNGDLQTSI) of human, rat and mouse SSTR2A Receptor coupled to KLH
Positive Control:	Human normal pancreatic islets
Synonyms:	anti-human, rat, mouse SSTR2 antibody, somatostatin receptor 2 (SSTR2, sst2, SSTR2A, sst2A) rabbit antibody raised against the carboxy-terminus of the human, rat and mouse SSTR2 receptor.
Localisation:	Cells and Surface Antigens

For reference only

For reference only
For Research Use Only. Not for Diagnostic or Therapeutic Use.