

Rabbit anti-Human Neurotensin Receptor 2 (NTR2) (AA 390-410) (CT)

#Cat: NB-49-006-50ul

Size: 50ul

Product Description:	NT2 antibody, human, C-term. NTR2 Neurotensin 2 receptor G Protein-coupled receptor Polyclonal antibody raised against the C-terminus of the human NT2 receptor
Immunogen:	residues 390-410 coupled to KLH (NP_036476) KLH- pqsptlmdtasgfgdppetrt
Host:	Rabbit
Form:	serum with 0.02% NaN ₃ as preservative
Quantity:	50 µl
Applications:	ICH (p), IC, W W = Western blot (1:5000) IC = Immunocytochemistry (1:5000) ICH (p) = Immunohistochemistry of formalin-fixed, paraffin- embedded tissues (1:2500)
Species cross-reactivity:	human
Sensitivity:	50 ng
Storage:	4 °C short term -20 °C long term



NT2 receptor immunoreactivity in a human gastric mucosa.