

## Bovine Plasma, K2-EDTA / Pooled, Female, Non-Sterile, Frozen

#Cat: NB-43-0001-10    Size: 10ml

#Cat: NB-43-0001-50    Size: 50ml

### Product Data Sheet

<b>Product</b>	Bovine Plasma, K2-EDTA / Pooled, Female, Non-Sterile, Frozen
<b>Product No</b>	NB-43-0001
<b>Lot</b>	31-112-042KW44/24
<b>Date</b>	30/10/2024
<b>Volume</b>	50 ml
<b>Presentation</b>	1 lot(s) of 50mL bovine plasma, frozen, pooled
<b>Species of Origin</b>	Cow
<b>Strain</b>	Holstein-Frisian
<b>Sex</b>	Female
<b>Remarks</b>	non medicated, non-immunized
<b>Anticoagulant</b>	K2EDTA
<b>Sterility</b>	aseptic
<b>Buffer</b>	n/a
<b>Preservative</b>	none
<b>Storage Conditions</b>	- 20°C.
<b>Expiry Date</b>	29/10/2025
<b>Disclaimer</b>	<p>This Material Should Be Handled as If Capable of Transmitting Infectious Agents. Please Use Universal Precautions.</p> <p>This Product Is Being Sold for Research And, Or Manufacturing Purposes Only. It Is Not to Be Used in Humans or Animals and Further Manufacturing That Will Result in A Finished Product That Contains Viable Leukocytes Is Prohibited. For In Vitro Use Only. The User Assumes All Responsibility for Its Usage and Disposal, In Accordance with All Regulations.</p>
<b>Country of Origin</b>	Germany