

## Rabbit anti-Human Neurokinin-1 Receptor (NK1R) (AA 387-407) (CT)

#Cat: NB-49-007-50ul

Size: 50ul

<b>Product Description:</b>	NK1 antibody, human, C-term. Substance P receptor SPR G Protein-coupled receptor Polyclonal antibody raised against the C-terminus of the human NK1 receptor
<b>Immunogen:</b>	residues 387-407 coupled to KLH (P25103) KLH- ssrsdsktmtesfsfnvls
<b>Host:</b>	Rabbit
<b>Form:</b>	serum with 0.02% NaN <sub>3</sub> as preservative
<b>Quantity:</b>	50 µl
<b>Applications:</b>	ICH (p), IC, W W = Western blot (1:5000) IC = Immunocytochemistry (1:5000) ICH (p) = Immunohistochemistry of formalin-fixed, paraffin- embedded tissues (1:2500)
<b>Species cross-reactivity:</b>	human
<b>Sensitivity:</b>	50 ng
<b>Storage:</b>	4 °C short term -20 °C long term



NK1 receptor immunoreactivity in a human prostate carcinoma.