
Rabbit anti-Human Neuromedin B Receptor (NMBR) (AA 342-361) (CT)

#Cat: NB-49-009-50ul

Size: 50ul

Product Description:	NMBR antibody, human, C-term. Neuromedin B receptor Neuromedin B-preferring bombesin receptor Bombesin receptor G Protein-coupled Receptor Polyclonal antibody raised against the carboxy-terminus of the human NMBR receptor.
Immunogen:	residues 342-361 coupled to KHL (CAH70473) KLH-nmvtnsvllnghsmkqemal
Host:	Rabbit
Form:	Antiserum
Quantity:	50 µl
Applications:	IC, W W = Western blot (1 µg/ml) IC = Immunocytochemistry (1 µg/ml)
Species cross-reactivity:	human
Storage:	4 °C short term -20 °C long term