

Goat Anti-Mouse IgG H+L antibody {AP}

#Cat: NB-22-113269-100

Size: 100uL

#Cat: NB-22-113269-200

Size: 200uL(2x100uL)

General Information

Product Type	Secondary antibodies
Short Description	
Applications	WB
Host/Source	Goat
Reactivity	Mouse

Product Properties

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	AP
Purification	Affinity purification
Dilution Range	WB:1:5000-1:50000
Formulation	PBS with 0.025% Sodium Azide, 0.75% BSA, 50% Glycerol, pH 7.3
Isotype	IgG
Storage Instruction	

Target Information

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Mouse IgG
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.