

Neo Horse Serum and Neo Donor Horse Serum

#Cat: NB-58-0005A	Size: 500ml
#Cat: NB-58-0005B	Size: 100ml
#Cat: NB-58-0006A	Size: 500ml
#Cat: NB-58-0006B	Size: 100ml

General Information

Neo Horse Serum and Neo Donor Horse Serum contain more immunoglobulins than Fetal Bovine Serum and have increased protein content. They are a cost-effective alternative to Fetal Bovine Serum.

Neo Donor Horse Serum features a controlled collection of the raw serum. In contrast to slaughtered horses, Neo donor horse serum has a lower endotoxin content and exhibits a lower level of batch to batch variation.

Features:

- Alternative to Calf Serum
- Suitable for ES cell culture
- Autologous cell culture system

Product	Volume	Cat. No.
Neo Horse Serum	100 ml	NB-58-0006B
	500 ml	NB-58-0006A
Neo Donor Horse Serum	100 ml	NB-58-0005B
	500 ml	NB-58-0005A

Product Specifications

Appearance	Clear amber liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened store at 4°C and use within 4-6 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.

Precautions and Disclaimer

This product is for research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (info@neo-biotech.com) or phone (+33 9 77 40 09 09).