

Anti-Mouse CD279/PD-1 antibody

Cat#NB-22-62351-50

Cat#NB-22-62351-100

Description

Rat monoclonal anti-Anti-Mouse CD279/PD-1 antibody is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunofluorescence

Product information

Product type	Primary antibodies
Applications	FC, WB, IHC, IF
Host/Source	Rat
Reactivity	Mouse
Clonality	Monoclonal
Clone ID	29F.1A12
Concentration	0.5 mg/mL

Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2a, k
Storage	Store at 2~8°C. Do not freeze.

Target

Gene ID	18566
Gene Symbol	Pdcd1
Uniprot ID	PDCD1_MOUSE

For reference only

For Research Use Only. Not for Diagnostic or Therapeutic Use.