Human Zonulin (Zonulin) Sandwich ELISA Kit

#Cat: NB-22-44016-24 Size: 24T
#Cat: NB-22-44016-48 Size: 48T
#Cat: NB-22-44016-96 Size: 96T

General Information

Product Type	Kits
Short Description	This Zonulin Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids.
Application	ELISA
Reactivity	Human


Product Properties

Detection	Colormetric
Sensitivity	0.47ng/mL
Detection Limit	0.78~50ng/mL
Assay Time	3.5h
Storage Instructions	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

Target Information

Sample Types	Serum, plasma and other biological fluids.
Gene Symbol	HP
Uniprot ID	HPT_HUMAN

For reference only



For reference only

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence.
This product must not be used as for diagnostic or other medical purposes.